

## Maxiprep library DNA

1. Thaw 2 vials of Stbl4 electrocompetent cells on ice
2. Spin down library DNA in tubes and transfer 2  $\mu$ l of each sub library to a new 1.5 ml tube
3. Add 54  $\mu$ l Stbl4 cells to each tube and incubate 10 minutes on ice
4. Transfer 25  $\mu$ l cells to a 1 mm gap cuvette and electroporate at EC1 (1.8 kV)
5. Immediately add 1 ml pre-warmed SOC and transfer to new 1.5 ml tube
6. Recover for 1-2 hour2 while shaking at 30 °C
  - a. Repeat transformation for the other half of the library
7. Take 10  $\mu$ l from each SOC tube and pool per sub set of the library in a new 1.5 ml tube
8. Add 80  $\mu$ l LB to each pool (total volume 100  $\mu$ l)
  - a. 10  $\mu$ l of this will give a  $10^{-3}$  dilution of sub library
9. Make 10-fold serial dilutions up to  $10^{-6}$
10. Plate 10  $\mu$ l of each dilution on a pre-warmed 10 cm ampicillin plate
11. Incubate all plate O/N at 37 °C
12. Transfer the remaining recovered cells to a 500 ml flask containing 100 ml amp medium and grow O/N while shaking at 30 °C
13. Maxiprep each sub library using the Qiagen Plasmid Plus Maxi Kit (according to the manufacturer's instructions) if the transformation efficiency exceeds the library size by 20-fold

## Generation library lentivirus

1. Seed 293T cells at  $7.5/10/12 \times 10^6$  per 15 cm plate and grow O/N at 37 °C in total 30.2 ml IMDM with 10% FCS without P/S
2. For transfection, choose the plate where confluency of the 293T cells is about 70%
3. Prepare transfection mix (amounts for two 15 cm plates):
  - a. 3 ml Optimem containing:
    - i. 13.4 µg pMD.G
    - ii. 26.8 µg pCMV dR8.91
    - iii. 20.2 µg Bassik whole genome library
  - b. 3 ml Optimem containing:
    - i. 181 µl TransIT293 transfection reagent
4. Mix solution a and b together
5. Incubate 30 min at RT
6. Add transfection mix dropwise to 293T cells
7. Remove supernatant from cells containing the virus 48 hours later
  - a. 293T cells should look very sick (membranes bubbly, etc) at this stage
8. Filter supernatant with 0.45 µm filter to remove 293T cells
9. Aliquot virus in 5 ml and 1.5 ml volumes
  - a. 5 ml for screen
  - b. 1.5 ml for titration of virus
10. Store virus at -80 °C

### Titration lentivirus for MOI = 0.3-0.5

For the screen each cell should only contain one sgRNA sequence. This experiment will determine the amount of virus to add to the target cell to get a multiplicity of infection (MOI) of 0.3-0.5.

1. Thaw 1.5 ml aliquot of virus at RT
2. Add  $2.5 \times 10^6$  HAP1 Cas9 to each well of a 6-well plate in total 5 ml complete IMDM
3. Add 0, 50, 75, 100, 150 or 200  $\mu$ l virus to a single 6-well
4. Spinfect:
  - a. Centrifuge plate 3,900 RPM for 1 hour at 37 °C
5. Transfer cells to incubator and grow for 48 hours at 37 °C
6. Wash cells with PBS
7. Add 500  $\mu$ l trypsin and incubate for 5-10 minutes at 37 °C
8. Add 4.5 ml complete IMDM
9. Resuspend cells completely and transfer 500  $\mu$ l to a new 6-well (2x)
10. To one well add 4.5 ml complete IMDM (-puromycin condition)
11. To the other well first add 2 ml complete IMDM and then 2.5 ml complete IMDM containing  $2 \times 3.5 \mu\text{g/ml}$  puromycin (+puromycin condition)
  - a. Repeat for all conditions
12. Grow cells for 48 hours at 37 °C
13. Wash cells with PBS
14. Add 500  $\mu$ l trypsin and incubate for 5-10 minutes at 37 °C
15. Add 4.5 ml complete IMDM
16. Count cells and determine MOI for each condition as follows:

$$\frac{\text{Cell count} + \text{puro}}{\text{Cell count} - \text{puro}}$$

$$\frac{\text{Cell count} + \text{puro}}{\text{Cell count} - \text{puro}}$$

- a. Choose virus amount for actual screen that will give you an MOI of 0.3-0.5

Note 1: when infecting the library for a screen make sure to exactly repeat this protocol to get the same MOI.

Note 2: with each new virus batch or library, repeat this process.